



FELINE THYROXINE ELISA TEST KIT

PRODUCT PROFILE AND INSTRUCTIONS

The T4 ELISA test is an immunoassay designed for the quantitative determination of thyroxine (T4) in serum/plasma samples of Feline and related species.

TEST PRINCIPLE

In the T4 ELISA Test, a certain amount of anti-T4 antibody is coated on microtiter wells. A measured amount of unknown sample and a constant amount of T4 conjugated with horseradish peroxidase are added to the microtiter wells. During incubation, the T4 and conjugated T4 compete for the limited binding sites on the anti-T4 antibody. After 60 minutes of incubation period, at 37C, the wells are washed 5 times with wash buffer to remove any unbound T4 conjugate. A solution of TMB is then added and incubated for 20 minutes, resulting in the development of a blue color. The color development is stopped with the addition of 2N HCl, and the absorbency is measured spectrophotometrically at 450nm. The intensity of the color formed is proportional to the amount of enzyme present and is inversely related to the amount of unlabeled T4 in the sample. A series of T4 standards assayed in the same way, a standard curve is constructed and the concentration of T4 in the unknown sample is quantified.

MATERIALS PROVIDED

1. Antibody-coated microtiter wells, 96-well plate
2. Reference Standard, 1 set, ready to use (0,5,10,25,50,100,300 ng/mL)
3. T4 HRP Conjugate Reagent, 12 mL
4. TMB Color Reagent (ready to use) 12 mL
5. Stop solution (2N HCl) 6mL
6. 20X Wash buffer, 20 mL
7. Instructions

MATERIALS REQUIRED, BUT NOT PROVIDED

1. Precision pipettes: 50uL, 100uL, 200uL, and 1.0mL
2. Disposable pipette tips
3. Vortex mixer or equivalent
4. Absorbent paper or paper towel
5. Graph paper
6. Microtiter plate reader

SPECIMEN COLLECTION AND PREPARATION

Serum should be prepared from a whole blood specimen obtained by acceptable medical techniques. This kit is for use with serum/plasma samples only.

STORAGE OF TEST KIT AND INSTRUMENTATION

Unopened test kits should be stored at 2-8^oC upon receipt and the microtiter plate should be kept in a sealed bag with desiccants to minimize exposure to damp air. Opened test kits will remain stable until the expiration date shown, provided it is stored as prescribed above. A microtiter plate reader with a bandwidth of 10nm or less and an optical density range of 0-2 OD or greater at a 450nm wavelength is acceptable for use in absorbency measurement.

REAGENT PREPARATION

1. All reagents should be brought to room temperature (25-28^oC) before use.
2. Must read and follow instructions carefully and any deviations may result in poor sensitivity or high background.
3. Prepare desired amount of wash buffer by diluting 1 part with 19 parts of distilled water. This buffer may be stored at 4-8C for 1-3 months.

ASSAY PROCEDURE

One must follow accurately these steps to ensure correct results. Use clean pipettes and sterile, disposable tips:

1. Secure the desired number of coated wells in the holder.
2. Dispense 25ul of standards, specimens, and controls into appropriate wells.

3. Dispense 100ul of T4 HRP-Conjugate Reagent into each well. Shake the plate for 30 seconds.. It is very important to shake the plate at this step.
4. Incubate at 37°C for 60 minutes.
5. Remove the incubation mixture by dumping plate contents into a waste container.
6. Rinse and dump the microtiter wells five (5) times with wash buffer.
7. Strike the wells sharply onto absorbent paper or paper towels to remove all residual water droplets.
8. Dispense 100ul of TMB color reagent into each well. Gently mix for 10 seconds.
9. Incubate at room temperature for 20 minutes, in the dark.
10. Stop reaction by adding 50ul of 2N HCl to each well.
11. Gently mix for 30 seconds. It is important to observe a color change from blue to yellow.
12. Read optical density at 450nm with a microtiter well reader.

Important note: The wash step is very critical and insufficient washing will result in poor precision and falsely elevated absorbency readings.

CALCULATION OF RESULTS

Calculate the mean absorbency value (A450) for each set of reference standards, specimens, controls and patient samples. Construct a standard curve by plotting the mean absorbency obtained from each reference standard against its concentration in ng/mL on graph paper, with absorbance values on the vertical or Y axis, and concentrations on the horizontal or X axis. Use the mean absorbency values for each specimen to determine the corresponding concentration of T4 in ng/mL from the standard curve.

EXPECTED VALUES AND SENSITIVITY

It is recommended that each laboratory should establish values to reflect differences specific to experimental conditions. The minimum detectable concentration of thyroxine by this assay is estimated to be 1 ng/mL.

REFERENCES

REFERENCES

1. Graves TK, Peterson ME. Diagnostic tests for feline hyperthyroidism Vet Clin North Am Small Anim Pract 1994 May;24(3):567-76
2. Brown RS, Keating P, Livingston PG, Bullock L. Thyroid growth immunoglobulins in feline hyperthyroidism. Thyroid. 1992 Summer;2(2):125-30.
3. Peterson ME, Kintzer PP, Cavanagh PG, Fox PR, Ferguson DC, Johnson GF, Becker DV. Feline hyperthyroidism: pretreatment clinical and laboratory evaluation of 131 cases. J Am Vet Med Assoc. 1983 Jul 1;183(1):103-10.
4. Peterson ME, Melian C, Nichols R Measurement of serum concentrations of free thyroxine, total thyroxine, and total triiodothyronine in cats with hyperthyroidism and cats with nonthyroidal disease J Am Vet Med Assoc 2001 Feb 15;218(4):529-36
5. Segal J, Troen BR, Ingbar SH. Influence of age and sex on the concentrations of thyroid hormone in serum in the rat. J Endocrinol 1982 May;93(2):177-81
6. Walker W.H.C. Introduction: An Approach to Immunoassay. Clin. Chem. 1977; 23: 384
7. Wisdom G.B. Enzyme-Immunoassay. Clin. Chem. 1976; 22: 1243
8. Hoffenberg R. Medicine 1978; 8: 392
9. Lieblich J., Utiger R.D. J. Clin. Invest. 1972; 51: 1939

Distributed by:

Xceltis

Xceltis GmbH
 Im Tal 12
 74909 Meckesheim / Germany
 Tel.: +49-6226-9724-18
 Fax: +49-6226-9724-19
 E-mail: info@xceltis.de

Revised 121002

FELINE T4 ELISA Test Kit

Professional use only

www.endocrinetech.com

35325 Fircrest Street, Newark, CA 94560-1003 * Phone (800) 745-0843 * (510) 745-0844 * Fax (510) 745-0977